

# TOKUYAMA® REBASE II

## Chairside Hard Denture Reline Material

Tokuyama® Rebase II is a chairside hard denture relining material. Without compromising any benefits from its predecessor, Tokuso Rebase, Tokuyama® Rebase II offers better polishability, lower heat generation for in-mouth curing, higher mechanical strengths, improved handling and less odor/taste. It cures in 4 minutes and 30 seconds intraorally. Tokuyama® Rebase II is indicated for hard denture relining and partial denture repairs.




### Benefits

- Lower chemical irritation
- Lower heat generation
- Durability
- Easier polishing
- One-appointment procedure
- Accurate adaptation of the denture
- Saves lab costs

### Indications

- Hard denture relining
- Partial denture repairs



# TOKUYAMA® REBASE II

Chairside Hard Denture Reline Material

## Product Packages Available

### Tokuyama® Rebase II Kit

- Tokuyama® Rebase II Powder/80g
- Tokuyama® Rebase II Liquid/50mL
- Adhesive/15mL
- Measuring cup (for Tokuyama® Rebase II Powder)
- Dropper (for Tokuyama® Rebase II Liquid)
- Rubber Cup (for mixing Tokuyama® Rebase II Powder and Liquid)
- Spatula (for spreading Tokuyama® Rebase II on the denture)
- Plastic Cup (for Adhesive)
- Brush (for Adhesive)
- Tokuso Resin Hardener II
- Spoon (for Tokuso Resin Hardner II)



### TOKUYAMA® REBASE II PRODUCTS

ITEM#	ITEM DESCRIPTION
20516	Tokuyama® Rebase II Kit
20556	Tokuyama® Rebase II Powder
20567	Tokuyama® Rebase II Liquid
20568	Adhesive
20172	Resin Hardner II
20133	Tokuso® Rebase Mr. Bond
23309	Sofreliner Tough® Kit

### Refill

- Tokuyama® Rebase II Powder/80g
- Tokuyama® Rebase II Liquid/50mL
- Adhesive/15mL
- Tokuso Resin Hardener II



### Tokuso® Rebase MR. Bond

MR. Bond is a primer for the bonding of self-curing acrylic resins with non-precious metals. It is a perfect companion to Tokuyama® Rebase II for metal plate denture relines.

### Sofreliner Tough®

Sofreliner Tough® is an addition-cured silicone chairside soft lining material for dentures. It has no odor, taste or chemical irritations. It offers a high tear resistance, high stain resistance and outstanding adhesion to the denture base. Sofreliner Tough® is designed to be used when durability is the primary concern. The paste cures in 5 minutes intraorally, and the kit comes with 2 different finishing points (shape adjusting point & finishing point).

